


## Clinical guide to Valgus ankle control

The most effective way of controlling ankle valgus is with a lateral strut design  
The ankle strap provides the 3rd point of support

1 medial wedging of the soft inset

Stages of support: 2 ankle strap

3 medial 1/2 SMO & ankle strap



## Recommended AFO's



### Posterior Lateral strut

foot drop  
PTTD  
recravatum  
ankle stability - valgus control



### Flexor Band AFO

foot drop  
for the active patient  
wide range of motion  
ankle stability - valgus control



### Lateral Strut

Foot drop  
PTTD  
ankle stability - valgus control



### Lateral Strut - Anterior Shell

provides knee stability at terminal stance  
ankle stability - valgus control  
foot drop