

flexorband afo

by CUSTOM COMPOSITE MFG, INC



The FlexorBand AFO is our first AFO with a joint. This lightweight, dorsiflexion assist AFO is designed for foot drop patients who need a wide range of motion at the ankle (-30 degrees dorsiflexion, unlimited plantarflexion).

The joint, located over the lateral malleolus, consists of integrated hemispherical stainless steel washers separated by a nylon washer. Tension is provided by a Gaffney-Flexor band anterior to the joint.

The FlexorBand AFO is custom made from a cast of the patient and comes with a foam liner, Velcro® strap and custom foam footbed.

For Casting Instructions, visit www.cc-mfg.com or call us toll free at 866-273-2230.

indications

Drop foot with wide ROM needed

recommended activity level

All activity levels, various sports (running, hiking, etc)

suggested I-codes

L-1980 Ankle foot orthosis, single, upright free plantar-dorsi flexion, solid stirrup, calf band/cuff, custom fabricated

L-2210 Addition to lower extremity, dorsiflexion assist (plantarflexion resist)

L-2755 Addition to lower-extremity orthosis, high strength, lightweight material, all hybrid lamination/prepreg composite, per segment, for custom-fabricated orthosis only

L-2820 Addition to lower extremity orthosis, soft interface for molded plastic, below knee section

